



## QuickTrace™ Mercury Analyzers Specifications of Standard Supplied Computer

Released 07-01-2010

Base Unit:	Dell OptiPlex 380 Minitower
Processor:	Pentium Dual Core E5400/ 2.7GHz, 2M 800FSB
Operating System:	Genuine Windows 7, XP PRO downgrade, English, Factory Install
Memory:	1 GB, Non-ECC, 1333MHz DDR2 1X1GB
Hard Drive:	160GB SATA, 3.0Gb/s Hard Drive with 8MB Data Burst Cache, Dell OptiPlex 780/580
DVD-ROM Drive:	16X DVD-ROM SATA, Data Only, Dell OptiPlex 780 Desktop or Minitower
Keyboard:	Dell USB Keyboard, No Hot Keys English, Black, OptiPlex
Mouse:	Dell USB 2-Button Optical Mouse with scroll, Black, OptiPlex
Video Card:	Integrated Video, GMA4500, Dell OptiPlex
Service:	<ul style="list-style-type: none"><li>✓ Dell Hardware Limited Warranty Plus Onsite Service Initial Year</li><li>✓ Dell Hardware Limited Warranty Plus Onsite Service Extended Years</li><li>✓ ProSupport for IT: 7X24 Technical Support for certified IT Staff, Initial Year</li><li>✓ ProSupport for IT: 7X24 Technical Support for certified IT Staff, 2 Year Extended</li><li>✓ ProSupport for IT: Next Day Parts and Labor Onsite Response Initial Year</li><li>✓ ProSupport for IT: Next Business Day Parts and Labor Onsite Response 2 Year Extended</li></ul>
Monitor:	Dell E190S, 19 Inch Flat Panel 19.0 Inch Viewable Image Size (3YR Limited Warranty)