

QuickTrace™ Mercury Analyzers

Specifications of Standard Supplied Computer

Released 06-01-2009

Base Unit:	OptiPlex 360 Minitower
Processor:	Pentium Dual Core E2200/ 2.2GHz, 1M 800FSB
Operating System:	Genuine Windows Vista Business Bonus - Windows XP PRO downgrade, English, Factory Install
Memory:	1 GB, Non-ECC, 667MHz DDR2 2X512MB
Hard Drive:	80GB SATA, 7200 RPM Hard Drive with Data Burst Cache, Dell OptiPlex 745 and 330
DVD-ROM Drive:	16X DVD-ROM SATA, Dell OptiPlex 360 Desktop or Minitower
Keyboard:	Dell USB Keyboard, No Hot Keys English, Black, OptiPlex
Mouse:	Dell USB 2-Button Optical Mouse with scroll, Black, OptiPlex
Floppy Disk Drive:	3.5 inch, 1.44MB, Floppy Drive Dell OptiPlex 740, 755, and 330
Video Card:	Integrated Video, GMA3100, Dell OptiPlex
Service:	<ul style="list-style-type: none">✓ Dell Hardware Limited Warranty Plus Onsite Service Initial Year✓ Dell Hardware Limited Warranty Plus Onsite Service Extended Years✓ ProSupport for IT: 7X24 Technical Support for certified IT Staff, Initial Year✓ ProSupport for IT: 7X24 Technical Support for certified IT Staff, 2 Year Extended✓ ProSupport for IT: Next Day Parts and Labor Onsite Response Initial Year✓ ProSupport for IT: Next Business Day Parts and Labor Onsite Response 2 Year Extended
Monitor:	Dell E178FP, 17 Inch Flat Panel 17.0 Inch Viewable Image Size OptiPlex, Precision and Latitude (3YR Limited Warranty)