



QuickTrace™ Mercury Analyzers Specifications of Standard Supplied Computer

Base Unit:	Dell OptiPlex 790 Minitower
Processor:	Pentium Dual CORE i3 2120 Processor (3.3GHz, 3M)
Operating System:	Microsoft Windows 7 Professional SP1, 32-Bit, English, Factory Install
Memory:	2GB,Non-ECC,1333MHz DDR3,1x2GB
Hard Drive:	250GB SATA 6.0Gb/s and 8MB Data Burst Cache
DVD-ROM Drive:	16X DVD-ROM SATA, Data Only
Expansion Slots:	1 PCI, 1 PCIe x16, 1 PCIe x16 (wired x4), 1 PCIe x1, 4 SATA
I/O Ports:	10 USB, 1 Serial, 1 LAN (integrated), 1 VGA, 1 DisplayPort
Keyboard:	Dell USB Keyboard, No Hot Keys English, Black
Mouse:	Dell MS111 USB 2-Button Optical Mouse with scroll, Black
Video Card:	Integrated Video,HD Graphics 2000
Service:	<ul style="list-style-type: none">✓ Dell Hardware Limited Warranty Plus Onsite Service Initial Year✓ Dell Hardware Limited Warranty Plus Onsite Service Extended Years✓ ProSupport for IT: 7X24 Technical Support for certified IT Staff, Initial Year✓ ProSupport for IT: 7X24 Technical Support for certified IT Staff, 2 Year Extended✓ ProSupport for IT: Next Day Parts and Labor Onsite Response Initial Year✓ ProSupport for IT: Next Business Day Parts and Labor Onsite Response 2 Year Extended
Monitor:	Dell E190S, 19 Inch Flat Panel 19.0 Inch Viewable Image Size (3YR Limited Warranty)